

Incubator Monitoring Solution


Affordable and simple laboratory and life science environmental monitoring with affordable sensors to monitor more critical conditions.

The Sonicu incubator monitoring sensor plugs seamlessly into the already robust monitoring system that includes temperature and air pressure differential monitoring.

The incubator monitoring solution includes temperature, humidity and CO₂, making it the ideal, cost-effective platform for life science professionals.

FEATURES

- Plug-and-Play Installation
- Simple mobile app configuration and monitoring
- Phone call alerts + text and e-mail alarms
- U.S.-based phone support
- Lifetime warranty with software subscription



MONITORING SPECIFICATIONS

CO₂ Accuracy/ Range	0 - 25% CO ₂ ; +/- (0.5% + 3% MV)
Temp Range/Accuracy:	-40 - 125°C; +/- 0.3°C
Humidity Range/Accuracy:	0 - 100% RH; +/- 2%



Sonicu incubator monitoring sensor