



# HEALTHCARE

The **Smart Choice** for ALL your monitoring needs.

**Sonicu provides superior operational visibility and uniform implementation across all departments and locations with one cost-effective solution. From a single vaccine refrigerator to an entire hospital campus, Sonicu's cloud-based platform provides a monitoring solution for every healthcare application.**



Temperature/cold storage and vaccine monitoring for medication/vaccine refrigerators, laboratory and research, cryogenic and ultralow, blanket warmers and food service.



Ambient temperature and relative humidity monitoring for operating suites, anesthetizing locations, sterile supply storage and sterilization areas.



Room conditions monitoring for critical- and non-critical environments; pressure, air velocity to calculate Air Changes per Hour (ACH).



Sound monitoring for improved patient safety and experience in NICUs, ICUs and patient areas.



Custom monitoring applications via Modbus/4-20 mA for connections between existing machines or a wide array of commercially available sensors.

Indiana University  
Health

BAYLOR  
HEALTH CARE SYSTEM  
FOUNDATION




LOMA LINDA UNIVERSITY  
HEALTH

Sonicu's Enterprise-wide solution consolidates all monitoring applications, including off-site locations, into one platform.


Sonicu's cloud-based platform enables real-time access to and viewing of all monitored points across an entire health system via a client-specific, customizable dashboard. Monitoring, alarming, logging and reporting are fully automated and standardized across your entire enterprise.

Sonicu wireless monitoring meets or exceeds compliance standards of The Joint Commission, HFAB, CDC, FDA, BOP, and VFC program.



# Why Sonicu?

## Safety, compliance, efficiency & savings.



### Safety:

Smart alarming sets Sonicu apart when it comes to protecting assets and the public. SMART alarm software features include:

- Custom thresholds with user-set high/low buffers to match monitored asset/equipment parameters for reduced alarm fatigue.
- Off-line alarm features such as batch alarming notification and auto-clear for improved operational efficiency.
- Alarm reminders to ensure all prior alarms have been addressed and cleared.
- Tiered, escalatable alarming to ensure proper notification goes to the right people at the right time.



### Effortless regulatory compliance:

- Consistent automated logging and auto-formatted reporting downloads streamline regulatory compliance.
- Intuitive software allows annotation and corrective action history to be included in regulatory reporting.
- Sonicu makes calibration plug-and-play. SNAP Calibration eliminates the downtime, expense and administrative worries of annual or recurring probe calibration.

### Save time and money:

- Recover lost hours and labor costs consumed by manual monitoring and hand logging.
- SaaS-based platform requires no IT resources, network equipment, additional software or maintenance costs.
- Plug-and-play capability virtually eliminates installation costs.

## All Sonicu monitoring solutions are:

Easily implemented and fully scalable without hardwiring or expensive retrofits.

Completely automated to eliminate manual monitoring, logging and reporting.

Totally secure and operate independently of existing IT infrastructure.

Warranted and supported by Sonicu's dedicated support team.

Cost-effective.

# Monitoring What Matters:

- Healthcare
- Pharmaceuticals
- Life Science
- Food Safety

"We're always looking for ways to improve our patient experience. Sonicu's technologies will help us do this and enable us to monitor temperatures, sound, humidity, and other environmental factors more efficiently."

Steve Long, CEO,  
Hancock Health  
and Hancock  
Regional Hospital.

[sonicu.com](http://sonicu.com)  
[info@sonicu.com](mailto:info@sonicu.com)

