

KN95 Masks

Sonicu KN95, 5-Layer PPE Masks

KN95 Masks

- $\geq 95\%$ Bacterial Filtration Efficiency
- 5 Protective Layers
- Adjustable Nose Bridge
- Elastic Ear Straps
- Soft and Comfortable
- Ergonomic Design



Part No. SMM-KN95
Operation Standard: GB2626-2006
Available in boxes of 25, or cases of 1,000

FDA Registered Product



Reusable mask protects against airborne particles including bacteria and viruses.

Non-Woven Fabric (50gm)

Melt-Blown Fabric (25gm)

Hot Air Cotton Fabric (50gm)

Melt-Blown Fabric (25gm)

Non-Woven Fabric (25gm)

Expiration Date: 2 years from manufacturing date.
Store at temperatures between -20°C and 38°C . Humidity $<80\%$.