

SoniCloud Monitoring Platform



Dashboard Access 24/7

Scaleable monitoring for any organization

Sonicu's cloud-based platform aggregates all monitored data in real-time to a single dashboard, with 24/7 access to live device status, plus historical data and reporting features from anywhere. Providing superior operational visibility enterprise-wide.



Cost-Effective, Wireless Monitoring Solutions

Temperature Monitoring



Sonicu's temperature monitoring protect vaccines, medications, research samples, reagents, food and more, 24/7.

Temperature Monitoring for:

- Hospitals, Clinics, Pharmacies
- Medical and Scientific Labs
- Cold Storage, IVR, Blood Banks
- Research Facilities, Cryogenic Storage
- Food Manufacturing, Distribution

Room Pressure/Condition



Room condition monitoring keeps critical spaces such as hospital rooms, pharmacies, operating rooms, and cleanrooms safe.

Pressure Monitoring for:

- Hospitals, Surgery Suites, Sterile Storage
- Clean Rooms, Compounding Pharmacies
- Scientific and Medical Research Labs
- Manufacturing and Industrial Facilities
- Data Centers, HVAC Systems

Sound Monitoring



Sound monitoring provides visual cues that prompt noise reducing behavior, and reports aggregate data to pinpoint noise.

Sound Monitoring for:


- NICU's for Neonatal Safety
- Patient Care Areas / Patient Satisfaction
- Sound Abatement for Construction
- OSHA Compliance for Industry

Full Visibility – Superior Asset Protection – Effortless Compliance

Sonicu is the leading provider for monitoring temperature, room pressure, sound, and other environmental conditions in the healthcare, pharmaceutical, food, life science, and other industries.

SoniCloud's SMART Tools Automate Regulatory Compliance Reporting.

Automated Reporting



Regulatory reports are automatically created including any annotations and corrective actions taken.

SMART Alarming

Users respond to alarms online. Records user, timeframe, asset type, application, and resolution action.

"SNAP" Calibration

Optional SNAP program ships NIST calibrated sensors to you just in time to replace expiring sensors.

Sonicu has the lowest cost-of-ownership, and provides the highest ROI.

- *Unlimited Users Included*
- *Software Updates Included*
- *Live Tech support Included*
- *Unlimited Data Storage Included*
- *Firmware Updates Included*
- *Unlimited Training Included*

Sonicu Plans and Included Features	Standard	Premium	Enterprise
Unlimited Data Storage (on the Cloud)	●	●	●
Unlimited Number of Users Allowed	●	●	●
Unlimited User Access to SoniCloud Academy	●	●	●
Unlimited Admin User Access to Sonicu Training	●	●	●
SMART Alarm Handler (with escalating alerts)	●	●	●
SMART Reporting Tools (auto-generated reports)	●	●	●
SoniCloud Hosting on AWS (Amazon Web Services)	●	●	●
Flexible Networking (WiFi, 4G/5G LTE, and Ethernet)	●	●	●
Dedicated Customer Success Manager	●	●	●
Live Customer Support (during business hours)	●	●	●
Hardware Warranty (with active subscription)	●	●	●
Diagnostic Information (signal strength / battery life)	●	●	●
Diagnostic Alarming*	○	●	●
Data Transmission at intervals of 5 min or greater	●	●	●
Data Transmission at any interval	○	●	●
Advanced Reporting (MKT, Excursion, etc.)	○	●	●
API Programming Access	○	●	●
API Programming Support	○	○	●
SSO Integration Support	○	○	●
Advanced SLA	○	○	●

* Under the Standard Plan, Diagnostic Alarming requires adding 1 point of monitoring to your subscription.

REV: 11-2020