

Healthcare Monitoring

Consistency and Ease Across
Your Enterprise Every Day.



Reliable Monitoring with Built-In Data Protection

From surgery suites to blood banks, Sonicu's monitoring solutions help your health organization stay ahead of everyday challenges while better protecting your patients, your reputation, and your bottom line.

System-Wide Visibility with Role-Based Access

Sonicu's monitoring solutions centralize environmental data onto **one platform**, giving managers a quick overview of all locations. Staff members see only what they need to with customized views based on their permission level.

Painless Calibration & Compliance

Sonicu sensors automatically upload **electronic calibration certificates** into SoniCloud, so you always have them on hand. Our SNAP Calibration Program takes away the burden of recalibration, and ships newly calibrated probes to your door.

Reliable Data Logging & Reporting

Recording data by hand is not only an inefficient use of labor hours, but can be a safety risk to patients. By automating the chore of data tracking and building out reports, your teams gain more bandwidth to focus on patient care, and your data records become **more consistent, accurate, and reliable**.

Real-Time Alerts on Your Phone

Alarm fatigue is a national safety hazard. Regulators (including Joint Commission) prioritize strategic alarm management as part of patient safety protocols. Sonicu's software, SoniCloud, **reduces alarm fatigue** using smart time buffers and tiered escalation, so you don't have to constantly verify readings. Sonicu will let you know if there's a problem.

Maximum, Redundant Data Protection

Sonicu's system provides multiple layers of redundancy to deliver **3X the data protection** and maintain a complete audit trail. Battery backup keeps Sonicu monitoring during power outages. Cellular failover prevents loss of connectivity if your primary network goes down. If connectivity is lost, our DataSync feature stores and forwards missing readings once you're back online.

Future-Proof Your Enterprise

Regulations change constantly, and violations can be costly. The key is to stay prepared. Facilities that upgrade to Sonicu gain access to a modernized platform, built for the regulatory demands of today and tomorrow.



Audit-Ready Reports

- Automated report delivery to teams via email—no manual exports
- On-demand or scheduled reports keep facilities audit-ready at all times
- Complete alarm history and corrective actions documented in every report

Complete Traceability


- Detailed audit trails capture and store time-stamped records of system activities to safeguard data integrity
- Digital documentation of corrective actions and staff involved to ensure team accountability
- Electronic calibration certificates upload automatically to the platform

Smooth Calibration Process

- NIST and ISO-Certified sensors that major compliance regulators approve and trust
- Digital calibration certificates upload automatically to the platform and are assigned to the matching meter
- All hardware information—meter IDs, sensor serial numbers, etc.—are uploaded to SoniCloud for you

Hassle-Free Re-Calibration

- Save on annual calibration costs when you enroll in Sonicu's SNAP Re-Calibration Program
- Newly calibrated probes delivered to your door when its time to replace your expiring probes
- No more manual calibrating and configuring—simply plug in your Sonicu sensor and SoniCloud will auto-update all the details to the matching meter



**View,
download,
& print
calibration
certificates
via SoniCloud**

Beyond Data Monitoring

Only Sonicu offers a complete monitoring and compliance solution on a single platform. Beyond monitoring, Sonicu delivers the reliability, connectivity, and control you need to manage operations across every department, room, and location.



Modern Mobile App

- **Instant alarms** sent via text, email, phone call, or push notification
- **Alert Management** allows teams to snooze, leave comments, or resolve alarms through the SoniCloud Mobile App
- **Tiered Escalation** minimizes alarm fatigue while enforcing team accountability and faster resolution

Unmatched Data Protection

- **Battery-Backup:** If the power goes out, Sonicu will continue monitoring, thanks to the built-in battery
- **Cellular Failover (optional):** During a network outage, Sonicu switches to connect to your preferred cell network for uninterrupted connectivity
- **DataSync:** In the event of a full outage, Sonicu will store all data until connectivity is restored, and will send the data to SoniCloud

Works With Your Tech Stack

- **Third Party Connectivity** pulls local equipment sensor data into SoniCloud
- **API Integration** streamlines Sonicu into existing healthcare workflow software
- **Modbus Adapter** available to integrate Particle Counting into the same platform

Why Choose Sonicu

Flexible Connectivity

Other providers expect you to connect their way. Sonicu supports Wi-Fi, Cellular, or Ethernet connectivity, so you can choose what fits your environment best, with built-in redundancy to keep you connected even during outages.



Intuitive & Ease of Implementation

Customize your SoniCloud site layout how it makes sense for your operations. Configure alarm thresholds, buffer times to prevent alarm fatigue, and tiered escalation to alert the right people at the right time.

Audit-Ready Reporting

Our reporting suite comes with pre-built reports, identifies trends and risk areas, and boosts accountability, so you can stay ahead of costly equipment failure without impacting daily operations.



Hendricks
Regional Health

"Our people were doing a solid job of manually logging temperatures but they also get called in a lot of directions and sometimes a log gets missed. And that's simply not acceptable because you always knew, in the back of your mind, that's a risk for the patients."

Peace of Mind is Just a Call Away



**Get back 2,500
Hours per Year**



**Save \$80k per
100 Hospital Beds**



**Maintain
99.99% Uptime**

Sonicu consolidates compliance documentation, audit reports, alarm management, and all hospital data monitoring under the same platform.

Our system combines wireless sensors, automated data logging, and built-in redundancy to deliver uninterrupted oversight of your facilities 24/7.

With over a decade of experience working with healthcare facilities, we have the knowledge and hardware necessary to best protect sensitive inventory and critical areas all while keeping facilities compliant with Joint Commission, FDA, and other healthcare regulations.

- **Your Network, Your Way:** Choose Cellular, Wi-Fi, or the SoniLink Hub for secure, reliable connectivity. Cell failover also available for backup during outages.
- **Intuitive Dashboards:** Our customizable views align with your workflow, with real-time insights accessible from any web-based device.
- **Compliance Made Easy:** With both automated and on-demand report generation built into our software, you'll always be ready for your next inspection.
- **Intelligent Alert Routing:** Our SMART Alarm Handler holds teams accountable with multiple alert methods, tiered escalation, and repeating alarm features.
- **Scalable System:** Our platform adapts as your business needs change without missing a beat.

**Set up a call with a Sonicu expert to
discuss your healthcare monitoring needs.**

Top Healthcare Facilities Choose Sonicu



OPTUM™



Piedmont
HEALTHCARE

vitalcare®
INFUSION SERVICES



OhioHealth



Indiana University Health



Community
Health Center
of Central Missouri

"When we put the numbers in front of our leadership and showed them how much time we were going to save with a remote system and automated monitoring versus the time we used to spend managing the system, it was a no-brainer."

Partnerships with Leading National GPOs and Healthcare Solutions Providers



PREMIER

Simplifi 797®

Wolters Kluwer



STERIS

vizient®

