


CO2 Incubator Monitoring Solution



Affordable and simple laboratory and life science environmental monitoring with affordable sensors to monitor more critical conditions.

The Sonicu incubator monitoring sensor plugs seamlessly into the already robust monitoring system that includes temperature and air pressure differential monitoring.

The incubator monitoring solution includes temperature, humidity and CO2, making it the ideal, cost-effective platform for life science professionals.

FEATURES

- Simple mobile app configuration and monitoring
- Phone call alerts + text and e-mail alarms
- U.S.-based phone support
- Lifetime warranty with software subscription

MONITORING SPECIFICATIONS

CO2 Range/ Accuracy	0 - 25% CO2 ; +/- (0.5% + 3% MV)
Temp Range/Accuracy	-40 - 125°C; +/- 0.3°C
Humidity Range/Accuracy	0 - 100% RH; +/- 2%
Dimensions	38.1 x 25.654 mm (1.5 in x 1.01 in)



Sonicu incubator monitoring sensor

