

One Platform. Open Integration. Smarter Workflows.

Sonicu's, open integration-ready platform securely connects monitoring data with your existing ecosystem. This gives teams a single source of truth while supporting the systems and workflows they already use.

Integrate Sonicu with reporting tools, analytics platforms, enterprise applications, or custom software for a complete, cohesive system.

Open API for Seamless Integration

Sonicu's open platform provides secure access to environmental monitoring data, alarms, device information, and system events, making it easy to integrate Sonicu however it makes sense for your operations.



Benefits of Sonicu's Open Platform

- **Connect critical data** across monitoring, reporting, analytics, and enterprise systems.
- **Automate operational workflows** by turning alarms and system events into timely actions.
- **Extend existing technology investments** without requiring a complete system replacement.
- **Standardize across the enterprise** with a secure, scalable platform designed for multiple facilities and applications.

Environmental Monitoring Data



Temperature



Humidity



Differential
Pressure



Water Leak



Power



Contact
Status

Sonicu Cloud Platform

Automated Monitoring and Compliance

Remote Alarming and Data Access

One Dashboard


API

Analytics Platforms Reporting Dashboards Enterprise Applications Custom Software

Technical Integration

Sonicu offers flexible integration options designed to fit your organization's architecture and reporting requirements. Whether you're building custom applications, integrating with enterprise software, or powering analytics dashboards, Sonicu provides secure, reliable access to your environmental monitoring data.

Access Key Platform Data



Integration Options

Secure REST API

Connect Sonicu with enterprise applications, custom software, and automation platforms using secure API access.

- Real-time access to platform data
- Secure authentication
- JSON responses
- Ideal for application integrations and custom development

SQL Read Replica

Feed Sonicu data directly into your existing system, granting you quick access to advanced analytics and custom reporting.

- Direct database access for reporting and analytics
- Optimized for business intelligence tools
- Supports enterprise reporting and data warehousing

Workflow Automation

Turn monitoring events into coordinated action by connecting Sonicu with the applications your teams use every day.

- Create service or maintenance tickets
- Route alarms to communication and escalation platforms
- Update enterprise dashboards
- Trigger compliance documentation
- Send data to analytics, quality, or facilities-management systems

Developer Resources

Complete documentation, authentication guides, and developer tools to help your team build integrations faster.

- **API Documentation:**
sonicumonitoring.com/api/documentation
- **Authentication & API Tokens:**
sonicumonitoring.com/api/pat